

ROTO SILENTA

BLOOD BANK CENTRIFUGE

FOR **CLINICAL & RESEARCH** APPLICATIONS

Hettich

CENTRIFUGES

The Hettich ROTO SILENTA floor-standing centrifuge is the incomparable choice in 12 Liter centrifuges on the market today.

Hettich VersaRotor™

Spin multiple bags, tubes and bottles with only one rotor, saving you time and money.

No windshield necessary

The ROTO SILENTA is more efficient and more user friendly simply because it does not need a windshield. Because of its powerful motor and dual sensor temperature probes, the SILENTA distributes its force and cooling more evenly around the sample itself.

Ergonomic

Without the need of a windshield the ROTO SILENTA is easier to load, easier to unload, easier to clean and ultimately easier on your back. The front facing control panel allows for easy entry of program settings, no need to lean over your unit to reach for controls.

Run-Silent™ Operation

The ROTO SILENTA has a decibel rating of ≤ 62 dB (as quiet as normal conversation).

Easy Maintenance

Without a wind shield getting in the way, the SILENTA is easy to clean and maintain by simply lifting the cover and wiping the surgical grade stainless steel bowl.

GMP/ GLP Compliant

Available with internal or external cooling for GMP and clean-room environments.

Compact

Surprisingly compact for its size, the ROTO SILENTA will take up less space in your lab than its competitors.

Installation and Support

Our commitment is to you. We are Here-to-HelpSM you with all your needs whether it be FREE installation and set-up or FREE applications support and training.



Find out more by visiting:

www.hettweb.com/roto-silenta

GSA Advantage![®]

GSA Contract Number: GS-07F-0586X

HETTICH ROTO SILENTA CENTRIFUGE


ADDITIONAL FEATURES

- Available with internal or external cooling for GMP and clean-room environments
- Non-wind shielded rotor for best sample temperature control
- Smooth operation with up to a 120 gram imbalance
- Control panel lock with key
- Free installation and training
- Bidirectional serial interface (optional)
- Bar code reader (optional)
- Hettinfo software package (optional)
- Run-Silent™ Operation (≤ 62 dB - as quiet as normal conversation)
- Five (5) year warranty


DESIGNED FOR:

Whole Blood
 Hard Spin
 Soft Spin
 Packed Red Blood Cells (PRBC)
 Platelet Rich Plasma (PRP)
 Platelet Concentrates
 Fresh/ Fresh Frozen Plasma (FFP)
 Cryoprecipitate/ Cryosupernat


BLOOD BAG PACKAGE - 6

	Tube Qty	Volume (mL)	RPM	G-FORCE	PACKAGE CONTENTS
	6	500	4,500	6,384	(1) Centrifuge (1) 90° Swing-out Rotor (6) Buckets (6) Adapters
	Pkg No.:		SILENTABLOODBAG6		
	Add Bio-containment Lids:		(-BC)		


BLOOD BAG PACKAGE - 6

	Tube Qty	Volume (mL)	RPM	G-FORCE	PACKAGE CONTENTS
	12	500	4,500	6,498	(1) Centrifuge (1) 90° Swing-out Rotor (6) Buckets (6) Adapters
	Pkg No.:		SILENTABLOODBAG12		

BOTTLE PACKAGE - 1L

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
	6	1,000	96 x 176	4,500	6,384	(1) Centrifuge (1) 90° Swing-out Rotor (6) Buckets
	Pkg No.:		SILENTABOTTLE1L			
	Add Bio-containment Lids:		(-BC)			

BOTTLE PACKAGE - 2L

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
	6	2,000	180 x 150 x 100	4,500	3,848	(1) Centrifuge (1) 90° Swing-out Rotor (6) Buckets
	Pkg No.:		SILENTABOTTLE2L			

The above packages are the most common packages available. More options and configurations are available. Please inquire.



GSA Contract Number: GS-07F-0586X

DISTRIBUTED BY:



HETTICH NORTH AMERICA
 Toll Free: 1-866-370-4388
 Email: info@hettweb.com

→ Discover more solutions online at: www.hettweb.com/centrifuges

Copyright © 2012 Hettich Instruments, LP. All rights reserved. For the most up to date information on Hettich please visit hettweb.com.
 Printed in the USA using soy based inks on FSC certified recycled paper (25% pre-consumer waste & 25% post consumer waste).

PS-ROTOSILENTA-E RevA 11-2012